



# VAC-TAK™

## THE VAC-TAK COMPOSITE INSTRUMENT

The Vac-Tak is a veneer placement instrument that makes it easy to handle and cement delicate porcelain or composite veneers. There are two sizes of suction tips – large, for maxillary teeth and small for mandibular teeth. There are also tips to fit the 3 mm light guide of a Microlux transilluminator.

- 1 autoclavable handle
  - Dropper water bottle
  - 15 large vacuum tips
  - 15 small vacuum tips
  - 15 large vacuum tip to fit 3 mm Microlux 2 light guide.
  - 15 small vacuum tips to fit 3 mm Microlux 2 light guide.
- The handle and vacuum tips are autoclavable.

Patent pending

43 Miry Brook Road • Danbury, CT 06810 • (203) 778-0200

**AdDent**  
INCORPORATED

Advancing Dental  
Excellence

ISO 9001 CERTIFIED & 13485 COMPLIANT | ALL PRODUCTS ARE FDA CLEARED

Made in the USA | [www.addent.com](http://www.addent.com)

## VAC-TAK VENEER PLACEMENT KIT


The Vac-Tak is a veneer placement instrument that makes it easy to handle and cement delicate porcelain or composite veneers. There are two size of suction tips – large, for maxillary teeth and small for mandibular teeth. It also has 2 sizes of vacuum tips to fit a 3mm light guide for our Microlux transilluminator.

### FEATURES


- Autoclavable handle
- Two size suction tips
- Autoclavable tips
- Silicone tips are transparent
- 3 mm lightguide tips

### BENEFITS


- Easy to hold ergonomic handle
- Firmly holds upper and lower veneers
- Silicone tips are reusable
- Tack cure through tips while holding veneer under pressure
- Hold and tack cure veneer with one instrument



Large silicone tip for upper teeth using Vac-Tak handle



Small silicone tip for lower teeth using Vac-Tak handle



Large silicone tip for upper teeth using 3 mm light guide



Small silicone tip for lower teeth using 3 mm light guide



PN: 660303 REV A 040124

Follow us on social media: [f](#) [in](#) [t](#) [v](#) [info@addent.com](mailto:info@addent.com)



**AdDent**  
INCORPORATED

Advancing Dental  
Excellence

ISO 9001 CERTIFIED & 13485 COMPLIANT | ALL PRODUCTS ARE FDA CLEARED



Made in the USA | [www.addent.com](http://www.addent.com)